

# HOUSEHOLD & INDUSTRIAL CARE



## Phosphonates

**ATMP & Salts** Amino-tris (methyleneephosphonic acid)

**ATMP n-Oxide** Amino-tris (methyleneephosphonic acid) n-Oxide

**DTPMPA & Salts** Diethylenetri-aminepenta (methyleneephosphonic acid)

**HEDP & Salts**-Hydroxyethylidene -1,1-diphosphonic acid

**PBTC & Salts** 2-Phosphonobutane -1,2,4,-tricarboxylic acid

## Others

1,2,3 Benzotriazole

4,4'-Bis-(2-sulfostyryl)-biphenyl disodium salt (BSBDS)

Alkylpolyglucosides

Glycolic acid 70%

2-Phenoxyethanol

## Chelating Agents

**EDTA 2Na** Ethylenediamine-tetraacetic acid disodium salt

**H-Quest**  
Sodium glucoheptonate

**MGDA 40 %** Methylglycine N,N-diacetic acid trisodium salt

**GLDA 38%** Tetrasodium N,N -bis(carboxylatomethyl) -L-glutamate

## Technical Preservatives

Various formulations available based on actives such as

**Bronopol**


**BIT** Benzisothiazolin-3-one

**DBNPA** Dibromonitrilopropionamide

**CIT/MIT** Isothiazolinones

**SDIC** Sodium dichloroisocyanurate

For your nearest service please visit: [www.connectchemicals.com](http://www.connectchemicals.com) · [inquiries@connectchemicals.com](mailto:inquiries@connectchemicals.com)



Phosphono carboxylate mix  
Tolytriazole  
Tributoxyethyl phosphate  
Sodium cumenesulfonate  
Silicone Fluids  
Surfactants